

NON-SURGICAL RHINOPLASTY :

Concise Method to Re-shape Nose Tips -Using REVOLAX™

Whether it's correcting a droopy tip or smoothing out a dorsal hump, non-surgical rhinoplasty (NSR) creates some of the most magical transformations in dermal filler aesthetics. Re-shaping the tip, however, does come with complications regarding the risk of tissue necrosis - an already long carried health warning of NSR. A less serious complication on the other hand is the creation of an aesthetically unappealing nasal tip - a commonly made mistake.

So, how do we reliably avoid both the risks of vascular injury and missing the aesthetic target?

Choosing a great product is key!


REVOLAX Deep and Sub-Q can both be used to create the perfect natural, non-surgical rhinoplasty. The thicker, long-lasting consistency of these products is perfect for augmentation.

Hitting your aesthetic target at the tip:

Injectors often struggle with precise bolus placement when augmenting the tip. You should be placing the bolus superficially, just below the dermis.

If you inject deep, it is less precise. It increases the likelihood of a boxy/bulbous tip and the potential to displace the underlying cartilage. The lower paired cartilage is held less firmly than the upper paired cartilage.

Remember that the upper nose is made of bone, the middle and lower parts of paired cartilage (upper and lower) with the septum



running down the middle like a rudder. Where the bone meets the upper paired cartilage is named the 'keystone' area. The upper and lower paired cartilage meet at the 'scroll' area. Deep injections above the scroll area do not tend to displace cartilage at the tip. However, deep injections do cause displacement of the lower cartilage (usually at the middle or medial crura), which is only correctable by dissolving the filler with hyaluronidase. This therefore should be avoided.

The scroll area is where the upper and lower paired cartilage meet. This is also approximately in line with the supra-tip breakpoint, or top of the alar groove.

Conclusion:

By staying superficial, mid-line and with a thorough understanding of the nose anatomy and areas to avoid - you will reshape nasal tips safely, successfully and with style.

Regarding 'how and when' to use a cannula, augmentation of the base, supra-tip and radix plus full access to our master tips - shop **Fox Academy's Advanced NSR E-learning** @ www.foxacademy.co.uk

Shop **REVOLAX** products from the UK licence holder, **Fox Group** @ www.foxpharma.co.uk | [@foxpharma_uk](https://twitter.com/foxpharma_uk)

